

PM-COMFORT™

PATIENTMATTRESS

IDÉ-RELAX-CARE

- Highcomfort-mattress
- Optimised for X-ray
- Hygienic and durable



PM-Comfort™ is designed for X-ray use. The mattress includes handles and lowfriction underside.



PM-Comfort™ is designed with unique pressure-relieving qualities. PM-Comfort™ can be delivered in different shape and size.

IDÉ-RELAX CARE is an Eson Comfort-catalogue of trademarks for hospitaluse!
UTGÅVA 05-06

PRODUCTINFORMATION

The mattress consists of a 15mm thick pressure-relieving upper layer. This layer is viscoelastic with a volume weight of ca 52kg/m³. Below this layer, there is a 15mm:s core with a volume weight of 55kg/m³. The mattress is flameresistant, in accordance with Swedish standards. The cover is designed to protect the mattress from urine and alcohol exposure. No CFC is used during production. PM-Comfort™ comes with a 1 year guaranty of function.

Standard versions

ART.NO	NAME
AICOMF	Angiomattress Width x Length x 5cm (No handles)
CTCOMF	CT-mattress 40 x 200 x 3cm (No handles)
PMCOMF	Patientmattress Width x Length x 3cm (No handles)
PMCOMFH	Patientmattress Width x Length x 3cm With handles

NOTE! Drawing or pattern, is always required for special shape or size.

Extras/spareparts

Special shape or size, handles, frictiontape to keep the mattress in place.

Care

The cover is cleaned with soap-alcohol solution or Virkon 1%.

CE Registrationcode 344-2757-35331

ESON COMFORT



ergonomiprofisen

ESON COMFORT

NO 1 in Patientcomfort

Centrumvägen 7 Box 105
SE-314 04 Landeryd

E-mail: info@esoncomfort.com
www.esoncomfort.com

Copyright 2003 Eson Comfort / Idé-Möbler i Landeryd AB. All rights reserved.

DISTRIBUTOR: